

## RifeMedic 5 Installation and connection:

### Computer requirements:

Laptop or Desktop, running on Windows XP or later and with an Audio speaker output.

### Installation, and test procedure:

#### 1. Install the software on your computer.

Observe during installation, the password appears in the License agreement text window.

Run the installed RifeMedic 5 program and follow the 3 steps to receive a pin number from me.

#### 2. When you have access to the program, test the software as follows:

Run the software and go to the Control page and click "Start" on Channel 1.

You should now hear the signal on your speakers. Click "Stop".

Do the same for Channel 2.

#### 3. Connect the audio cable from your computer audio output to the back of the RifeMedic 5.

Connect the power cable.

Turn the two level knobs fully anti-clockwise to minimum.

On your computer, adjust the speakers volume to maximum.

Switch the RifeMedic on. The LCD will display channel 1 and channel 2 intensity and load.

On the software, go to the Programs page.

Double-click on any program name to load it in the list.

For channel 1, on the left hand, on the bottom of the window, click:

Clear/Stop, then Load Selected, then Goto Control.

On the control page, click the Start button for channel 1.

On the RifeMedic, adjust channel 1's intensity knob clockwise and observe the 'Output Signal' leds.

The output leds will increase in brightness as you increase the output intensity.

You will normally have the hand electrodes connected to the RifeMedic, and while holding them in both hands, adjust the intensity to a level where the sensation of the signal is comfortable or where you may just feel the signal.

Stop the running program by clicking the 'Stop' button on your computer screen.

Adjust the RifeMedic intensity knob anti-clockwise, back to minimum output intensity.

Repeat the exercise for Channel 2.

Please read the training session for details on all the adjustments and features.

We offer a comprehensive bio-resonance / Rife frequency therapy course from time to time to interested persons.

Best regards,

Dr. Pieter van Wyk

011 794 2750 078 800 5753 modeltec@icon.co.za

[www.rifemedic.co.za](http://www.rifemedic.co.za) [www.effata.co.za](http://www.effata.co.za) [www.innovationinstruments.co.za](http://www.innovationinstruments.co.za)